

Lifeguard

Guarding Lives Innovatively...



Lifeguard-BPH-01

NIBP Holter

Lifeguard

Features:

- Transfer data to PC-based Software via USB port, accurate analysis systems.
- 2.4" TFT display, 320x240 pixel.
- Up to 300 items of memory for common user and 358 for Ambulatory Blood Pressure Monitoring user
- Clinically proven accuracy
- Automatic cuff inflation and deflation
- Reliable power supply with 2 AA batteries
- Battery status indicator
- Portable design
- Measurement Range: 0mmHg-280mmHg
- Accuracy: ± 3 mmHg
- Resolution: 1 mmHg
- Measurement item: SYS/DIA/MAP/PR

Technical Specifications:

- Measurement method: Oscillometry
- Operation model: Manual/Automatic/STAT
- Measurement unit: mmHg/KPa selectable
- Measurement types: Systolic, Diastolic, Mean
- Measurement range:
 - Adult pressure: 10 to 270 mmHg; Pulse rate: 30 to 220;
 - Pediatric pressure: 10 to 200 mmHg; Pulse rate: 30 to 220;
 - Neonate pressure: 10 to 135 mmHg; Pulse rate: 30 to 220
- Resolution: 1 mmHg
- Over-pressure protection: Yes
- Input/output: USB2.0
- Accuracy: Pulse rate: $\pm 2\%$ or ± 3 BPM which is greater
- Blood pressure : Meets ANSI/AAMI SP10-1992
- Safety: Independent safety microprocessor monitors cuff pressure, measurement time as well as the operation of the main microprocessor

Software Screen shorts:



Lifeguard

Medical Systems Pvt. Limited

An ISO 9001:2008 Certified Co.

SCO 83, 1st Floor, Sector 40-C, Chandigarh, INDIA
Tel. : +91-172-260016 Fax: +91-172-4606016 Mobile: +91-76962 57100
mail : sales@lifeguard-med.com

www.lifeguard-med.com