

Lifeguard

Lifeguard
Guarding Lives Innovatively...



Lifeguard-O2
Lifeguard-O2 Sleep Monitor



Features :

- Continuous monitoring of SpO2 and Pulse rate.
- 10 hours continuous recording.
- Monitoring on mobile device via Bluetooth.
- Overnight monitor of SpO2 and Pulse rate.
- Vibration is triggered by low oxygen.
- Screening for OSA (OBSTRUCTIVE SLEEP APNEA).

Specifications:

Weight	35g (main unit)	Size	44x25x15 mm (main unit)
Display	LED	Wireless	Bluetooth 4.0 BLE
SpO2 range	70% to 100%	Pulse rate	30 to250 bpm
Battery	Rechargeable Lithium-polymer Battery		
Charge Time	Less than 2 hours to 90%		
Vibrator	Trigged by low oxygenation		
Recording Time	Sleep Mode: 10 hours		
Recorded parameters	SpO2, Heart Rate, Steps, Motion		

Lifeguard

Medical Systems Pvt. Limited

An ISO 9001:2008 Certified Co.

www.lifeguard-med.com

SCO 83, 1st Floor, Sector 40-C, Chandigarh, INDIA
Tel. : +91-172-260016 Fax: +91-172-4606016 Mobile: +91-76962 57100
mail : sales@lifeguard-med.com