

# Lifeguard

*Guarding Lives Innovatively...*

# Lifeguard



**Lifeguard-BP**  
Smallest Personal BP Device

## Technical Specification:

- Size 120mmx60mm(Dia)
- Weight 160g
- Pressure range 20-280mm Hg
- Pressure accuracy +/-3 mmHg or 2%, whichever is greater
- Pulse rate range 30 to 250 bpm
- PR accuracy +/-2bpm or +/-2%, whichever is greater
- Cuff Upper arm circumference 22-32cm
- Inflation Manual
- Water proof IP22
- Drop 1.5m
- Power USB Chargeable, one charge 500 tests
- Wireless Bluetooth 4.0 BLE
- APP Android / iOS
- Package Main unit, cuff, charge cable, manual

## **Lifeguard** **Medical Systems Pvt. Limited**

An ISO 9001:2008 Certified Co.

SCO 83, 1st Floor, Sector 40-C, Chandigarh, INDIA

Tel. : +91-172-260016 Fax: +91-172-4606016 Mobile: +91-76962 57100

mail : sales@lifeguard-med.com

[www.lifeguard-med.com](http://www.lifeguard-med.com)