

Lifeguard

Guarding Lives Innovatively...

Lifeguard



Lifeguard-ECG-03

3 Channel ECG Machine

Features :

- Memory of 680 exams
- Recording Paper size: 80mmx20m
- 7 Inch Color LCD for Waveform Preview
- ECG data transfer to PC for storage, editing and for common A4 paper printing.
- Auto & Manual Operation
- Filter Out Powerful anti-interference AC and HUM filter
- Liquid Crystal Display
- Built-in Rechargeable Battery
- Full Power Capacity ensures of at least 50 Patients' ECG Recording

Technical Specifications:

Lead	Standard 12 Lead	Rhythm Lead	Standard 12 Lead
Input Circuit	Isolated; Protection against Pace making & defibrillation	Operation Mode	AUTO/MANU
Input Impedance	≥50 MΩ	Frequency Response	0.05-200Hz(-3dB)
Sensitivity	2.5, 5, 10, 20mm/mV	Paper Speed	5,12.5, 25, 50mm/s±3%
Filter	AC: 50 Hz(-20dB) EMG: 35-45Hz(-3dB)	Recording Paper	Roll Paper 80mmx20m
Skin Voltage Tolerance	±600mV	Sampling Rate	1ms
CMRR	≥110dB	A/D Converter	12 bit
Baseline Control	Automatic Position	LCD	7 graph LCD
Printer	Thermal array printer	Safety Level	Class I, Type CF
Weight	4.3Kgs	Power Requirements	AC: 220 V, 50Hz/110V,60Hz DC: Built in rechargeable battery Supporting 50 exams after complete charging

Lifeguard

Medical Systems Pvt. Limited

An ISO 9001:2008 Certified Co.

SCO 83, 1st Floor, Sector 40-C, Chandigarh, INDIA

Tel. : +91-172-260016 Fax: +91-172-4606016 Mobile: +91-76962 57100

mail : sales@lifeguard-med.com

www.lifeguard-med.com