

OOPSIE HEROES

Oopsie Heroes is a fun & simple bedwetting alarm for kids. It combines an app with a tiny sensor that teaches children how to make it to the toilet on time at night. Goodbye Oopsie moments!



Features:

- The tiny sensor tracks those Oopsie moments, finds the pattern and helps your child learn to wake up on time.
- The sensor is super simple to use and can be easily fastened to all undergarments and pj's. No special clothing required.
- After an Oopsie moment, the sensor will record and send this information to the Oopsie Heroes App.
- Instead of Bluetooth, the Oopsie sensor uses a simple sound to send the information to the app.

OOPSIE HEROES APP

- The Oopsie Heroes app automatically logs Oopsie moments. This helps you as parents spot the pattern and understand when your child needs support, motivation or a high five after a dry night. This helps to motivate your child to become dry in record time.
- Included features are a rewarding system, Parent Voice Alarm and tracking of progress.

DRY IN RECORD TIME

- Oopsie Heroes will train the brain of your child to recognize the bladder signals. Over time the child will learn to recognize these signals and wake up on their own.
- Within 6-12 weeks of using Oopsie Heroes your child will be free from nighttime accidents. Dry nights, better days!